



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-CD274

RNASCOPE® PROBE INFORMATION


Catalog number:	600861 (Manual)*, 600866 (VS assay)*, 600867 (LS assay)* 600869 (2.5 VS assay)*, 600868 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	CD274 molecule (CD274)
Alternate names:	B7-H; B7H1; PDL1; PD-L1; PDCD1L1; PDCD1LG1
Entrez Gene ID:	29126
Accession#:	NM_014143.3
Unigene:	Hs.521989
Target Sequence:	Homo sapiens CD274 molecule (CD274), transcript variant 1, mRNA
Probe Region:	124 - 1122
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE cervical cancer tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-CD274 probe (Right; Cat. No. 600861), 40X



RNAScope® Manual Brown Assay performed on FFPE renal cell carcinoma tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-CD274 probe (Right; Cat. No. 600861), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/29126>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.

Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.