



TARGET PROBE: HS-CMET


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION

Catalog number:	310051(Manual)*, 310056 (VS assay)*, 310057 (LS assay)* 310059 (2.5 VS assay)*, 310058 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	MET
Alternate names:	HGFR; AUTS9; RCCP2; c-Met; DFN97
Entrez Gene ID:	4233
Accession#:	NM_001127500
Unigene:	Hs.132966
Target Sequence:	Homo sapiens met proto-oncogene (hepatocyte growth factor receptor) (MET), transcript variant 1, mRNA
Probe Region:	1236-2257
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown assay performed on FFPE colorectal cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-CMET probe (Bottom; Cat. No. 310051), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/4233>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.