


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION

Catalog number:	411251 (Manual)*, 411256 (VS assay)*, 411257(LS assay)* 411259 (2.5 VS assay)*, 411258 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	CX3CR1
Alternate names:	V28; CCRL1; GPR13; CMKDR1; GPRV28; CMKBRL1
Entrez Gene ID:	1524
Accession#:	NM_001337.3
Unigene:	Hs.78913
Target Sequence:	Homo sapiens chemokine (C-X3-C motif) receptor 1 (CX3CR1) transcript variant 4, mRNA
Probe Region:	1608-2554
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE liver cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-CX3CR1 probe (Bottom; Cat. No. 411251), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/1524>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.