


RNASCOPE® PROBE INFORMATION

Catalog number:	311321 (Manual)*, 311326 (VS assay)*, 311327 (LS assay)* 311329 (2.5 VS assay)*, 311328 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	CXCL13
Alternate names:	BLC; BCA1; ANGIE; BCA-1; BLR1L; ANGIE2; SCYB13
Entrez Gene ID:	10563
Accession#:	NM_006419.2
Unigene:	Hs.100431
Target Sequence:	Homo sapiens chemokine (C-X-C motif) ligand 13 (CXCL13), mRNA
Probe Region:	30-1186
ZZ pair number:	16

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE gastric cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-CXCL13 probe (Bottom; Cat. No. 311321), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/10563>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.