



TARGET PROBE: HS-FGFR2


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION

Catalog number:	311171 (Manual)*, 311176 (VS assay)*, 311177 (LS assay)* 311179 (2.5 VS assay)*, 311178 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	FGFR2
Alternate names:	BEK; JWS; BBDS; CEK3; CFD1; ECT1; KGFR; TK14; TK25; BFR-1; CD332; K-SAM
Entrez Gene ID:	2263
Accession#:	NM_000141.4
Unigene:	Hs.533683
Target Sequence:	Homo sapiens fibroblast growth factor receptor 2 (FGFR2), transcript variant 1, mRNA
Probe Region:	1743-3033
ZZ pair number:	25

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE lung cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-FGFR2 probe (Bottom; Cat. No. 311171), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/2263>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.