


RNASCOPE® PROBE INFORMATION

Catalog number:	310791 (Manual)*, 310796 (VS assay)*, 310797 (LS assay)* 310799 (2.5 VS assay)*, 310798 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	FGFR3
Alternate names:	ACH; CEK2; JTK4; CD333; HSFGR3EX
Entrez Gene ID:	2261
Accession#:	NM_000142.4
Unigene:	Hs.1420
Target Sequence:	Homo sapiens fibroblast growth factor receptor 3 (FGFR3), transcript variant 1, mRNA
Probe Region:	536-2798
ZZ pair number:	28

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE bladder cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-FGFR3 probe (Bottom; Cat. No. 310791), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/2261>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.