


RNASCOPE® PROBE INFORMATION


Catalog number:	311181 (Manual)*, 311186 (VS assay)*, 311187 (LS assay)* 311189 (2.5 VS assay)*, 311188 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	NRG1
Alternate names:	GGF; HGL; HRG; NDF; ARIA; GGF2; HRG1; HRGA; SMDF; MST131; MSTP131; NRG1-IT2
Entrez Gene ID:	3084
Accession#:	NM_013956.3
Unigene:	Hs.453951
Target Sequence:	Homo sapiens neuregulin 1 (NRG1), transcript variant HRG-beta1b, mRNA
Probe Region:	442-2977
ZZ pair number:	45

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE lung cancer tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-NRG1 probe (Right; Cat. No. 311181), 40X



RNAScope® Manual Brown Assay performed on FFPE colorectal cancer tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-NRG1 probe (Right; Cat. No. 311181), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/3084>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.