



TARGET PROBE: HS-PDGFRB


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION

Catalog number:	548991 (Manual)*, 548996 (VS assay)*, 548997 (LS assay)* 548999 (2.5 VS assay)*, 548998 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	PDGFRB
Alternate names:	IMF1; IBGC4; JTK12; PDGFR; CD140B; PDGFR1; PDGFR-1
Entrez Gene ID:	5159
Accession#:	NM_002609.3
Unigene:	Hs.509067
Target Sequence:	Homo sapiens platelet-derived growth factor receptor, beta polypeptide (PDGFRB), mRNA
Probe Region:	523-2984
ZZ pair number:	40

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE renal cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-PDGFRB probe (Bottom; Cat. No. 548991), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/5159>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.