



TARGET PROBE: HS-VTCN1


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION

Catalog number:	418081 (Manual)*, 418086 (VS assay)*, 418087 (LS assay)* 418089 (2.5 VS assay)*, 418088 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	V-set domain containing T cell activation inhibitor 1
Alternate names:	B7X; B7H4; B7S1; B7-H4; B7h.5; VCTN1; PRO1291
Entrez Gene ID:	79679
Accession#:	NM_001253849.1
Unigene:	Hs.546434
Target Sequence:	Homo sapiens V-set domain containing T cell activation inhibitor 1 (VTCN1) transcript variant 2 mRNA
Probe Region:	745 - 1844
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE ovarian cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-VTCN1 probe (Bottom; Cat. No. 418081), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/79679>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.