



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-IDO1

RNASCOPE® PROBE INFORMATION


Catalog number:	602681 (Manual)*, 602686 (VS assay)*, 602687 (LS assay)* 602689 (2.5 VS assay)*, 602688 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	IDO1
Alternate names:	IDO; INDO; IDO-1
Entrez Gene ID:	3620
Accession#:	NM_002164.5
Unigene:	Hs.840
Target Sequence:	Homo sapiens indoleamine 2,3-dioxygenase 1 (IDO1), mRNA
Probe Region:	24 - 1081
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE ovarian cancer tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-IDO1 probe (Right; Cat. No. 602681), 40X



RNAScope® Manual Brown Assay performed on FFPE esophageal cancer tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-IDO1 probe (Right; Cat. No. 602681), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/3620>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.