



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-PDCD1

RNASCOPE® PROBE INFORMATION

Catalog number:	602021 (Manual)*, 602026 (VS assay)*, 602027 (LS assay)* 602029 (2.5 VS assay)*, 602028 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	PDCD1
Alternate names:	PD1; PD-1; CD279; SLEB2; hPD-1; hPD-I; hSLE1
Entrez Gene ID:	5133
Accession#:	NM_005018.2
Unigene:	Hs.158297
Target Sequence:	Homo sapiens programmed cell death 1 (PDCD1), mRNA
Probe Region:	107 - 1238
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE ovarian cancer tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-PDCD1 probe (Right; Cat. No. 602021), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/5133>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.