



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-CD70

RNASCOPE® PROBE INFORMATION

Catalog number:	419331 (Manual)*, 419336 (VS assay)*, 419337 (LS assay)* 419339 (2.5 VS assay)*, 419338 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	CD70
Alternate names:	CD27L; CD27LG; TNFSF7
Entrez Gene ID:	970
Accession#:	NM_001252.4
Unigene:	Hs.501497
Target Sequence:	Homo sapiens CD70 molecule (CD70) mRNA
Probe Region:	10 - 885
ZZ pair number:	13

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE lymphoma tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-CD70 probe (Bottom; Cat. No. 419331), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/970>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.