



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-FOXP3

RNASCOPE® PROBE INFORMATION

Catalog number:	418471 (Manual)*, 418476 (VS assay)*, 418477 (LS assay)* 418479 (2.5 VS assay)*, 418478 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	FOXP3
Alternate names:	JM2; AIID; IPEX; PIDX; XPID; DIETER
Entrez Gene ID:	50943
Accession#:	NM_014009.3
Unigene:	Hs.247700
Target Sequence:	Homo sapiens forkhead box P3 (FOXP3) transcript variant 1 mRNA
Probe Region:	568 - 1620
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE gastric cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-FOXP3 probe (Bottom; Cat. No. 418471), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/50943>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.