



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-GNLY

RNASCOPE® PROBE INFORMATION

Catalog number:	407371 (Manual)*, 407376 (VS assay)*, 407377 (LS assay)* 407379 (2.5 VS assay)*, 407378 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	GNLY
Alternate names:	LAG2; NKG5; LAG-2; D2S69E; TLA519
Entrez Gene ID:	10578
Accession#:	NM_006433.3
Unigene:	Hs.105806
Target Sequence:	Homo sapiens granulysin (GNLY) transcript variant NKG5 mRNA
Probe Region:	4 - 716
ZZ pair number:	11

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE lymphoma tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-GNLY probe (Bottom; Cat. No. 407371), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/10578>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.