



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-PRF1

RNASCOPE® PROBE INFORMATION

Catalog number:	407381 (Manual)*, 407386 (VS assay)*, 407387 (LS assay)* 407389 (2.5 VS assay)*, 407388 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	PRF1
Alternate names:	P1; PFP; FLH2; PFN1; HPLH2
Entrez Gene ID:	5551
Accession#:	M28393.1
Unigene:	Hs.2200
Target Sequence:	Human perforin mRNA complete cds
Probe Region:	7 - 1587
ZZ pair number:	19

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE gastric cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-PRF1 probe (Bottom; Cat. No. 407381), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/5551>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.