



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-CCL5

RNASCOPE® PROBE INFORMATION

Catalog number:	549171 (Manual)*, 549176 (VS assay)*, 549177 (LS assay)* 549179 (2.5 VS assay)*, 549178 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	CCL5
Alternate names:	SISd; eoCP; SCYA5; RANTES; TCP228; D17S136E; SIS-delta
Entrez Gene ID:	6352
Accession#:	NM_002985.2
Unigene:	Hs.514821
Target Sequence:	Homo sapiens chemokine (C-C motif) ligand 5 (CCL5), mRNA
Probe Region:	47 - 1223
ZZ pair number:	16

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE Esophageal cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-CCL5 probe (Bottom; Cat. No. 549171), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/6352>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.