



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-CCL7

RNASCOPE® PROBE INFORMATION

Catalog number:	425261 (Manual)*, 425266 (VS assay)*, 425267 (LS assay)* 425269 (2.5 VS assay)*, 425268 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	CCL7
Alternate names:	FIC; MARC; MCP3; NC28; MCP-3; SCYA6; SCYA7
Entrez Gene ID:	6354
Accession#:	NM_006273.3
Unigene:	Hs.251526
Target Sequence:	Homo sapiens chemokine (C-C motif) ligand 7 (CCL7) mRNA
Probe Region:	9 - 769
ZZ pair number:	11

* Target region is the same across manual and automation

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE ovarian cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-CCL7 probe (Bottom; Cat. No. 425261), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/6354>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.