



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-IFNB1

RNASCOPE® PROBE INFORMATION

Catalog number:	417071 (Manual)*, 417076 (VS assay)*, 417077 (LS assay)* 417079 (2.5 VS assay)*, 417078 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	IFNB1
Alternate names:	IFB; IFF; IFNB; IFN-beta
Entrez Gene ID:	3456
Accession#:	NM_002176.2
Unigene:	Hs.93177
Target Sequence:	Homo sapiens interferon beta 1 fibroblast (IFNB1) mRNA
Probe Region:	2 - 795
ZZ pair number:	16

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE lymphoma tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-IFNB1 probe (Bottom; Cat. No. 417071), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/3456>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.