



TARGET PROBE: HS-IL4


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION

Catalog number:	315191 (Manual)*, 315196 (VS assay)*, 315197 (LS assay)* 315199 (2.5 VS assay)*, 315198 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	IL4
Alternate names:	BSF1; IL-4; BCGF1; BSF-1; BCGF-1
Entrez Gene ID:	3565
Accession#:	NM_000589.3
Unigene:	Hs.73917
Target Sequence:	Homo sapiens interleukin 4 (IL4), transcript variant 1, mRNA
Probe Region:	1-529
ZZ pair number:	16

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE pancreatic cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-IL4 probe (Bottom; Cat. No. 315191), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/3565>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.