



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-IL7

RNASCOPE® PROBE INFORMATION

Catalog number:	424251 (Manual)*, 424256 (VS assay)*, 424257 (LS assay)* 424259 (2.5 VS assay)*, 424258 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	IL7
Alternate names:	IL-7
Entrez Gene ID:	3574
Accession#:	NM_001199888.1
Unigene:	Hs.591873
Target Sequence:	Homo sapiens interleukin 7 (IL7) transcript variant 4 mRNA
Probe Region:	55 - 1482
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE esophageal cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-IL7 probe (Bottom; Cat. No. 4242511), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/3574>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.