



TARGET PROBE: HS-IL10


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION

Catalog number:	602051 (Manual)*, 602056 (VS assay)*, 602057 (LS assay)* 602059 (2.5 VS assay)*, 602058 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	IL10
Alternate names:	CSIF; TGIF; GVHDS; IL-10; IL10A
Entrez Gene ID:	3586
Accession#:	NM_000572.2
Unigene:	Hs.193717
Target Sequence:	Homo sapiens interleukin 10 (IL10), mRNA
Probe Region:	122 - 1163
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE ovarian cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-IL10 probe (Bottom; Cat. No. 602051), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/3586>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.