



TARGET PROBE: HS-IL12A


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION

Catalog number:	402061 (Manual)*, 402066 (VS assay)*, 402067 (LS assay)* 402069 (2.5 VS assay)*, 402068 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	IL12A
Alternate names:	P35; CLMF; NFSK; NKSF1; IL-12A
Entrez Gene ID:	3592
Accession#:	NM_000882.3
Unigene:	Hs.673
Target Sequence:	Homo sapiens interleukin 12A (natural killer cell stimulatory factor 1, cytotoxic lymphocyte maturation factor 1, p35) (IL12A), mRNA
Probe Region:	136 - 1297
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE renal cell carcinoma tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-IL12A probe (Bottom; Cat. No. 402061), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/3592>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker. Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.