



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-IL17A

RNASCOPE® PROBE INFORMATION

Catalog number:	310931 (Manual)*, 310936 (VS assay)*, 310937 (LS assay)* 310939 (2.5 VS assay)*, 310938 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	IL17A
Alternate names:	IL17; CTLA8; IL-17; CTLA-8; IL-17A
Entrez Gene ID:	3605
Accession#:	NM_002190.2
Unigene:	Hs.41724
Target Sequence:	Homo sapiens interleukin 17A (IL17A), mRNA
Probe Region:	23-1102
ZZ pair number:	19

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE gastric cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-IL17A probe (Bottom; Cat. No. 310931), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/3605>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.