



ADVANCED CELL DIAGNOSTICS


TARGET PROBE: HS-IL20

RNASCOPE® PROBE INFORMATION

Catalog number:	412201 (Manual)*, 412206 (VS assay)*, 412207 (LS assay)* 412209 (2.5 VS assay)*, 412208 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	IL20
Alternate names:	IL-20; IL10D; ZCYTO10
Entrez Gene ID:	50604
Accession#:	NM_018724.3
Unigene:	Hs.272373
Target Sequence:	Homo sapiens interleukin 20 (IL20) mRNA
Probe Region:	4 - 1229
ZZ pair number:	20

* Target region is the same across manual and automation

PROBE DATA



RNAscope® Manual Brown Assay performed on FFPE bladder cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-IL20 probe (Bottom; Cat. No. 412201), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/50604>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.
Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.