



TARGET PROBE: HS-KIT


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION

Catalog number:	606401 (Manual)*, 606409 (2.5 VS assay)*, 606408 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	KIT
Alternate names:	PBT; SCFR; C-Kit; CD117
Entrez Gene ID:	3815
Accession#:	NM_001093772.1
Unigene:	Hs.479754
Target Sequence:	Homo sapiens v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog (KIT), transcript variant 2, mRNA
Probe Region:	364-1325
ZZ pair number:	20

* Target region is the same across RNAScope® manual and automation assays

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE breast cancer tissue under standard condition, using the dapB Negative Control probe (Top; Cat. No. 310043) and Hs-KIT probe (Bottom; Cat. No. 606401), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/3815>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.

Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.