


ADVANCED CELL DIAGNOSTICS

RNASCOPE® PROBE INFORMATION


Catalog number:	451231 (Manual)*, 451239 (2.5 VS assay)*, 451238 (2.5 LS assay)*
Species:	Homo sapiens
Target Gene:	TUG1
Alternate names:	TI-227H; LINC00080; NCRNA00080; TUG1
Entrez Gene ID:	55000
Accession#:	NR_002323.2
Unigene:	Hs.554829
Target Sequence:	Homo sapiens taurine up-regulated 1 (non-protein coding) (TUG1) transcript variant 3 long non-coding RNA
Probe Region:	1028 - 2009
ZZ pair number:	20

* Target region is the same across RNAScope® manual and automation assays

PROBE DATA



RNAScope® Manual Brown Assay performed on FFPE gastric cancer tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-TUG1 probe (Right; Cat. No. 451231), 40X



RNAScope® Manual Brown Assay performed on FFPE colorectal cancer tissue under standard condition, using the dapB Negative Control probe (Left; Cat. No. 310043) and Hs-TUG1 probe (Right; Cat. No. 451231), 40X

For more information regarding this target gene: <http://www.ncbi.nlm.nih.gov/gene/55000>

Note: ACD always recommends RNA quality verification using control probes, prior to running a target marker.

Expression pattern may vary from data shown here based on tissue preparation methods and/or associated sample biology.